

Glossary

BED WIDTH AND LENGTH

Standard senior living beds are 36" wide. A wider 39", 42", or 48" bed provides more space to sleep and reposition, which increases comfort and could reduce falls. They may feature:

- A specific, fixed width and/or length
- Built-in adjustments allowing for multiple width options
- Detachable accessories that allow the bed to expand to other widths
- Adjustable lengths

HEIGHT/TRAVEL RANGE

Indicates both the low height and high height of the bed.

DECK STYLE

Good decks provide a firm sleeping surface that won't change over time. The most common styles:

- Slat**
Metal slats provide a strong, supportive frame and are easy to clean.
- Grid**
Large-diameter wire welded to a frame provides high mattress ventilation and eliminates sharp edges.
- Spring**
Usually consists of a thick wire grid attached to a metallic frame by springs.

Pan
Offers the best orthopedic support, has no sharp edges, is easy to clean, and comes in two styles:

- Perforated, which allows for high mattress ventilation
- Composite, which allows for ease in removal and cleaning

LOCKING MECHANISMS


Step Lock
An easy-to-use safety lock that can be initiated by the foot to provide immediate immobilization of the bed.

Floor Lock
Typically extended to place the rubber foot of the lock onto the floor surface to create drag or friction to immobilize the bed. The placement and number of locks is determined by the caster type.

Standard Caster Locks
These locks are individually applied to each caster.

AUTOMATIC CONTOUR

Activated by simultaneously pressing the head rise/lower and foot rise/lower buttons on the bed control to initiate a head and foot cradle effect. Some bed controls have a default automatic contour feature triggered by pressing one button.



CONTROL LOCKOUT

A resident safety feature that allows the caregiver to lock the hand pendant functionality.

MOBILITY

Enables bed mobility via the bed's wheel locking mechanism, depending on the high/low position of the bed. For example, some beds are only mobile when the bed is in its lowest position, while others are mobile at any height.

TRENDELENBURG POSITION

The Trendelenburg position is a patient position that is used for certain kinds of medical procedures. When a patient is placed in this position, he or she is lying face up, and the surface on which the patient lies is angled so that the feet are above the head. In the reverse position, the head lies above the feet.

See our full-line Catalog V65, pages 1762-1771, and our Healthcare Catalog V10, pages 130-139, for our complete selection of beds.



BRAND	PRODUCT NAME	WEIGHT CAPACITY	DECK LENGTH	DECK WIDTH	ADJUSTABLE LENGTH	TRAVEL RANGE	DECK STYLE	LOCKING MECHANISM	LOCKS AT ANY HEIGHT	MOBILITY AT ANY HEIGHT	INCLUDES HEAD/FOOT BOARDS	INCLUDES RAILS	AUTOMATIC CONTOUR	# OF MOTORS	CONTROL LOCKOUT AVAILABLE	RAIL TYPE	WARRANTY	PART #
 Yes, you can.®	CS9 FX600	600 Lb	80"	36, 39, Or 42"	No	9 To 30"	Slat	Step Lock	Yes	Yes	Yes	No	Yes	3	Yes	Assist Device, Assist Rail Or Thinksoft Device	2-Year On Electronics, 5-Year On Mechanical; Verify Serial Number	Search the following SKUs on hdsupplysolutions.com : 677420, 677421, 677422, 677423, 677425, 677426, 677427, 677428 Online Only
	CS7	500 Lb	76 To 80"	36, 39, Or 42"	Yes	7 To 30"	Slat	Floor Lock System	Yes	Yes	Yes	Optional	Yes	3	Yes	Assist Bar, Assist Rail, Or Separate Thinksoft Positioning Device	2-Year On Electronics, 5-Year On Mechanical, 10-Year On Frame; Verify Serial Number	See V10, Page 134, For Configured Options
	CS5	500 Lb	76 To 80"	36, 39, Or 42"	Yes	9 To 30"	Slat	Step Lock	Yes	Yes	Yes	Optional	Yes	4	Yes	Assist Bar, Assist Rail, Or Separate Thinksoft Positioning Device	2-Year On Electronics, 5-Year On Mechanical, 10-Year On Frame; Verify Serial Number	See V10, Page 135, For Configured Options
	CS3	450 Lb	76 To 80"	36"	Yes	8-1/2 To 27"	Slat	Step Lock	No	No, Only In Lowest Position	Yes	Optional	Yes	4	Yes	Assist Bar, Assist Rail, Or Separate Thinksoft Positioning Device	2-Year On Electronics, 5 Years On Mechanical, 10-Year On Frame; Verify Serial Number	See V10, Page 136, For Configured Options
	SC900DLX	450 Lb	76 To 80"	36"	Yes	9-7/8 To 20-1/4"	Grid	Standard Caster Locks	Yes	No, Only In Lowest Position	Yes	Yes	No	3	No	Assist Rail	10-Year On Frame And Welds, 2-Year On Mechanical, 1-Year On Electronics	673507, 673508: See V10, Page 138
	820DLX	450 Lb	76 To 80"	36"	Yes	13-1/2 To 27"	Grid	Standard Caster Locks	Yes	No, Only In Lowest Position	Yes	Yes	No	3	No	Assist Rail	10-Year On Frame And Welds, 2-Year On Mechanical, 1-Year On Electronics	673502, 673503: See V10, Page 138
	UltraCare	500 Lb	76 To 80"	35, 39, Or 42"	Yes	7 To 30"	Grid	Ultra Lock	Yes	Roll In Low Only Or Roll At Any Height Options	Yes	Optional	Yes	4	Yes	Assist Device Or Deluxe Assist Handle	Lifetime On Welds, 15-Year On Frame, 5-Year On Parts And Components	See V10, Pages 130-131, For Configured Options
	EasyCare	500 Lb	76 To 80"	35, 39, Or 42"	Yes	7 To 30"	Grid	Step Lock	Yes	Roll In Low Only Or Roll At Any Height Options	Yes	Optional	Yes	5	Yes	Assist Handle, Assist Device, Or Deluxe Assist Handle	Lifetime On Welds, 15-Year On Frame, 5-Year On Parts And Components	See V10, Page 132-133, For Configured Options
	Elite Riser	600 Lb	76 To 80"	36"	Yes	7 To 30"	Grid	Z-Lok Mechanism	Yes	Yes	No	No	Yes	3	Yes	Pivot Assist Bar, Half Length Assist Device, Or Fixed Assist Bar	15-Year On Frame, 5-Year On Motor And Control Box, 3-Year On Other Electronics	See V10, Page 137
	Maxi-Rest	750 Lb	80"	48"	No	15-1/2 To 29-12"	Pan	Standard Caster Locks	Yes	Yes	Yes	Yes	Yes	3	No	Head End Side Rails	15-Year On Frame, 5-Year On Motor And Control Box, 3-Year On Other Electronics	676428: See V10, Page 61 Or 139
	Bari Low EC Bed	850 Lb	84"	36, 42, Or 48"	No	9-3/4 To 26"	Pan	Standard Caster Locks	Yes	Yes	Yes	Yes	Yes	3	No	Head End, 1/4 Length	5-Year On Frame, 1-Year On Electrical Parts And Components	679296: See V10, Page 61 Or 139
	Sentinel Bariatric Low Bed	1000 Lb	86"	48"	No	9-1/4 To 30"	Pan	Individual Locking Casters	Yes	Yes	Yes	Yes	Yes	3	No	Head And Foot End Side Rails	5-Year On Frame, 1-Year On Electrical Parts And Components	676502: See V10 Pages 61 Or 139

View More Bed Configurations Online

Visit hdsupplysolutions.com and search **healthcare beds & mattresses** to see other bed configurations that meet your needs. Contact your sales representative if you have any questions.



Joerns EasyCare



NOA Elite Riser



Gendron Sentinel Bariatric Low Bed